

## Quantitative Rapid Test Detecting Aflatoxin M1 In Milk

AflaSensor is a rapid assay in dipstick format detecting quantitatively the contamination of milk samples by Aflatoxin M1 molecules.

The test can be used in-situ for daily controls or in laboratories for the analysis of series of milk samples. This test does not require any sample processing, cleaning or extraction.



### BENEFITS



**QUANTITATIVE**



**RAPID**

Results in 10 minutes.



**SENSITIVE**

Adapted to specific regulation (EU and US).



**RELIABLE AND ROBUST**



**USER-FRIENDLY**

Easily performed on site or in lab.



**COST-EFFECTIVE**

### RESULTS

The use of Readsensor 2 allows an accurate quantification of the Aflatoxins contamination.



### SENSITIVITY

TOXIN	EU MRL	US ACTION LEVEL (SL)
	KIT041 (QUANTITATIVE)	KIT078 (QUANTITATIVE)

Aflatoxin M1 (PPT)


Detection range from 30 to 150

Detection range from 200 to 750

## HOW TO USE




1  
Mix milk  
and reagents



2  
Place the vial  
to heating block



3  
Close and  
push "START"



4  
Place dipstick  
to Heatsensor

INCUBATE **10** MINUTES



5  
Interpret  
results

QUANTITATIVE

# AflaSensor Quanti

## PRODUCTS REFERENCES

### KIT041v2

AFLASENSOR QUANTI MK 0,05 PPB - 96 PC

### KIT078

AFLASENSOR MK 0,5 PPB 96 pc

APP032 HEATSENSOR DUO

APP008 AIRSHIELD DUO

APP052 HEATSENSOR OCTO

APP033 AIRSHIELD OCTO



APP088

READSENSOR 2



**unisensor**  
DIAGNOSTIC ENGINEERING

Liège Science Park  
Rue Louis Plescia, 8  
B-4102 Ougrée (Liège)  
BELGIUM

[www.unisensor.be](http://www.unisensor.be)

Phone +32 4 335 30 00  
[info@unisensor.be](mailto:info@unisensor.be)

Reseller