


Imagine how technology  
could transform your milk testing


**aurox**  
NEXT-GENERATION MILK TESTING

Embrace automation and  
set a new standard in milk testing



# Goodbye Pipettes! Milk Testing Made Simple...

1



Insert a reagent disc

2



Pour milk into the funnel

3



Start the analysis and  
get the results of your  
sample as fast as under  
3 minutes\*

4



Track and manage the  
results and related  
information



Discover the future of dairy safety with aurox® by Unisensor – the first fully automated system that combines lateral flow immunoassay with fluid handling technology. The aurox® eliminates the need for extensive training and manual pipetting, while enhancing testing robustness and repeatability.

The aurox® streamlines dairy operations by delivering real-time data and verified antibiotic test results whether during transport or in laboratory analysis.




#### Automated sample handling

- No more pipetting
- No microwells, no manual mixing
- No volume concerns
- Self cleaning



#### High test performance

- High throughput
- Accurate
- Reproducible



#### Result analysis

- Digital result interpretation
- User-friendly interface
- Unbiased, confident decisions



#### Integrated data management

- Seamless operation, zero downtime
- 4G, Wi-Fi, Ethernet connectivity with remote monitoring
- Unisensor data management or LIMS integration

**ILVO** Validation ensures compliance with the highest industry standards for reliability and performance.



Designed and Manufactured in Belgium.

The aurox® DISCB201 detects **beta-lactams (including cefalexin)** and **tetracyclines** in milk (cow, goat and sheep), meeting EU residual limits

- **Rapid Testing:** Get results for 1 sample in as fast as under **3 minutes\***
- **High Throughput:** Processes up to **6 tests simultaneously** in less than 5 minutes\*
- **Compact format:** Packaged in a disc format with up to **12 tests per disc**, for easy handling and compact storage

\*Timing based on tests with milk at room temperature; sample loading time not included.



# Unisensor - Who We Are

## Our Mission

Bringing smart diagnostic solutions from the laboratory to the field.

We have consistently met field challenges by innovating with advanced, intelligent diagnostic solutions that go beyond traditional laboratory boundaries.

With our unmatched technologies, we don't simply strive for product enhancement, we aspire to revolutionize the contours of our industry.

## Our Vision: Peace Of Mind in Food

Consumer confidence is essential for food producers. With our products, we help them guarantee superior quality food that meets consumers' expectations.



Liège Science Park • Rue Louis Plescia, 8 • 4102 Ougrée (Liège) • Belgium

Phone +32 4 252 66 02 • [info@unisensor.be](mailto:info@unisensor.be) • [www.unisensor.be](http://www.unisensor.be)